

DVO-1000MD

Medical DVD Recorder



Overview

The simplicity and power of DVD recording - an essential part of your medical imaging system

Sony now makes all the advantages of the DVD format available to clinicians who need to record high-quality images. Fully compliant with latest Medical Devices Directive (MDD) standards, the DVO-1000MD is a DVD Recorder that has been designed for use in a wide range of healthcare environments. It allows recording and playback of images from a wide range of visible and non-visible light modalities. Compact, ruggedly-engineered and easy to use, the DVO-1000MD offers all the benefits of removable disc media over S-VHS tape. These include digital image quality and near-instant browsing and access to files, plus low media costs and reduced library space requirements.

This equipment is intended for use by qualified medical professionals only.

Exceptional image clarity

The DVO-1000MD uses industry-standard MPEG-2 video compression. This allows images to be recorded with high clarity that's required in medical applications while optimising the use of disk space.

Quick random access by title and chapter

You can quickly access a desired title or chapter simply by selecting it from a list that is automatically created and instantly displayed simply by pressing the 'LIST' button. These powerful random access capabilities speed up the process of reviewing recordings after and examination.

Instant recording and extra security

Thanks to the built-in 80GB Hard Disk Drive (HDD), recording begins the moment the REC button is pressed. What's more, the HDD provides even more security for your recordings because you won't lose data in the unexpected event of power outage in recording or damage onto the DVD disc.

Low cost and wide interoperability of media

Adopting the inexpensive DVD+RW format of removable disc media, the DVO-1000MD offers playback compatibility with a wide range of devices. The DVD+RW media is low-cost and eliminates the need for a finalisation process after recording. It is also ideal for long-term archiving, as well as giving obvious space-saving benefits compared to VHS tapes.

Patient data handling

The DVO-1000MD facilitates the recording of patient data on the disc, allowing simpler data management.

Features

Long recording time

The DVO-1000MD provides a long recording time of up to 3 hours on a single disc. Three recording modes are available, depending on the picture quality and

recording duration required :- HQ Mode : 60 minutes SP Mode : 120 minutes LP Mode : 180 minutes

Quick random access by title and chapter

The DVO-1000MD makes records as 'Titles' - recording segments from when the REC button is pressed until the STOP button is pressed. Within a Title, you can set 'Indexes' to images of interest by pressing the PAUSE button during recording. A segment between two index points is handled as a 'Chapter'. When playing back the disc, quick access to the top of a title or chapter is possible with the press of a button.

High speed colour picture search

The DVO-1000MD offers +/-3 times normal speed playback as standard. Using an optional control panel such as SVRM-100A or DSRM-10, +/-10 times normal speed playback becomes available with easy-to-use Jog/Shuttle dial operation. Furthermore, a maximum +/-30 times normal speed playback is possible from other medial equipment connected using the RS-232C protocol.

Variety of remote controls

For remote control from other medical equipment, the DVO-1000MD supports the RS-232C protocol - widely used in the Sony SVO-9500MD S-VHS Medical VTR. Three remote control units are also available as options : SVRM-100A, DSRM-10 and FS-20. A USB 2.0 (Full Speed) interface is also provided.

Wide range of inputs

The DVO-1000MD is equipped with a wide range of interfaces including analogue composite, S-Video and analogue audio I/Os and an audio monitor output. An i.LINK (DV) input is also provided, facilitating high-quality digital video recording with minimal signal degradation.

PAL/NTSC switchable

The DVO-1000MD can be operated for both PAL (625) and NTSC (525) systems by a simple menu setting.

Specifications

General

Power requirements	AC 100V to 240V, 50/60 Hz
Dimensions	212 x 128.5 x 382 mm
Mass	6kg
Recording Format	Video: MPEG-2 compression Audio: Dolby Digital format
Recording Media	DVD+RW (2.4x and 4x speed)
Medical Standards	EN60601-1, DN60601-1-2, UL60601-1, CSA C22.2 No.60601.1

Input/Output

Analogue Composite Input	BNCx2, with loop-through
S-Video Input	4-pin DINx2, with loop-through
Analogue Composite Output	BNCx1
S-Video Output	4-pin DIN x 1

Input/Output

Analogue Audio Input	RCA Pin x 2 (L/R)
Analogue Audio Output	RCA Pin x 2 (L/R)
Monitor Audio Output	Monitor RCA Pin x 1
i.LINK (DV IN)	i.LINK 6-pin x 1, IEEE1394
Remote Control	RS-232C x 1, D-sub 9-pin
	Remote 1 x 1, stereo mini jack (for SVRM-100A and DSRM-10 controllers)
	Remote 2 x 1, stereo mini jack (for FS-20 foot switch)
	USB 2.0 x 1 (full speed)

Supplied Accessories

Operation Manual

Protocol Manual

DPW120MD

DVD+RW Blank Disc

Canopus DVD InfoEditor

Editing Software

Gallery



© 2004 - 2026 Sony Corporation. All rights reserved. Reproduction in whole or in part without written permission is prohibited. Features and specifications are subject to change without notice. The values for mass and dimension are approximate. All trademarks are the property of their respective owners.