

## LMD-KT10

Calibration Kit for LMD-DM series  
Medical Diagnostic Displays



### Overview

This easy-to-use PC software and sensor kit supports consistent display quality control with Sony LMD-DM series monitors – an essential requirement with medical diagnostic display applications.

Compatible with LMD-DM50, LMD-DM30, LMD-DM20, LMD-DM30C, LMD-DM20C medical diagnostic displays.

### Features

#### **Display Calibrator**

Allows easy, accurate display calibration with supplied calibration sensor.

#### **Allows easy, accurate display calibration with supplied calibration sensor.**

Allows accurate measurement of display luminance.

#### **Display Quality Controller**

Provides acceptance and consistency tests according to AAPM, DIN and JESRA quality assurance guidelines.

#### **Display Utility Software**

Allows central control and management of display information by remote system administrator.

## Gallery

