

## UPC-24SA

Self-laminating colour printing pack (S-size) for UP-25MD and UP-D25MD medical printers



### Overview

A6 (S-size) self-laminating dye sublimation colour printing pack for use with Sony UP-25MD and UP-D25MD digital colour printers.

Pack contains 180 sheets of print media and ribbon (60 sheets x 3 packs).

Media size 100 x 90 mm.

### Specifications

#### Specifications

Type of Media	Self-Laminating Color Printing Pack
Contents	180 sheets of print paper 3 roll of ink ribbon
Paper/Film Size	100 x 90 mm (A7) 4 x 3 1/2 inches
Print Pages for One	180 sheets of print paper (60sheets x 3)

Roll

3 roll of ink ribbon

---

Supported Printers

UP-25MD, UP-D25MD

---

## Gallery

