

## UPC-510

Colour printing pack for UP-50, UP-D50 and UP-51MD medical printers.



### Présentation

A5 dye sublimation colour printing pack for use with Sony UP-50, UP-D50 and UP-51MD printers.

Pack contains 200 sheets of print media plus ink ribbon.

Media size 148 x 210 mm.

### Spécifications techniques

#### Spécifications techniques

Type of Media	Color Printing Pack
Contenu	200 sheets of print paper 1 roll of ink ribbon
Paper/Film Size	210mm x 148 mm (A5) 8 3/8 x 5 7/8 inches
Print Pages for One Roll	200 sheets of print paper 1 roll of ink ribbon
Supported Printers	UP-51MD,UP-51MDP,UP-D50

Informations quant à la loi européenne sur les données : [Cliquez ici](#) pour vérifier si votre produit et/ou service associé Sony est soumis à la loi européenne sur les données.

### Galerie



modifiées sans préavis. Les valeurs pour le poids et les dimensions sont approximatives. Les autres marques appartiennent à leurs propriétaires respectifs.