

CU-C31

Condenser cardioid Mic Capsule
for use with DWM-02 and DWZ
series handheld transmitters



Overview

Designed for vocal applications including live music performance, the CU-C31 is a condenser cardioid mic capsule designed for use with the DWM-02 and DWZ series handheld transmitters

Condenser Microphone Capsule for use with DWM-02

Handheld Transmitter

Specifications

General Specifications

Capsule Type	Condenser
Dimensions	1 7/8" Diameter x 2 1/2" Height
Directivity	Cardioid
Frequency Response	60 Hz - 20 kHz
Weight	4.8 oz.

Gallery

