

## HVO-AUDIOKIT

Audio narration kit for the  
HVO500/550MD Medical Grade  
DVD Recorder



### Overview

#### **Audio narration kit for the HVO500/550MD Medical Grade DVD Recorder**

This kit will give the HVO500/550MD an audio feature enabling an audio track to be recorded onto a DVD as video is being recorded. This feature can later be used for review, and training purposes.

#### **Highlights**

- Rane Microphone Amplifier
- Sony ECM-673 Microphone
- XLR/RCS Adapter
- XLR male to XLR female cable
- RCA/RCA Y Cable

## Gallery

