

HVO-KEYBOARD

USB Keyboard for HVO-500/500MD



Overview

USB Keyboard for HVO-500/500MD

This HVO-KEYBOARD is designed specifically for the medical market. It features a completely sealed water-resistant robust rubber casing. The keyboard can easily be washed down to prevent the potential spread of harmful bacteria. Patient data can be entered efficiently by utilizing a standard QWERTY keyboard layout with tactile feedback keys.

Highlights

- Sealed Robust Rubber Keyboard
- USB cable is included
- Full Transport control of all front panel controls of HVO-500/550MD

Gallery

