

## PWSK-4501

Internal Memory Array (2TB)



### Overview

#### **Internal Memory Array (2TB)**

The PWSK-4501 Internal Memory Array (2TB) is for the PWS-4500 live server which is equipped with 2TB storage capacity at standard configuration. You can expand storage capacity up to 8 TB by installing three PWSK-4501.

#### **Expanding the storage capacity of the PWS-4500 up to 8 TB.**

By installing three boards of PWSK-4501 in the PWS-4500, a large storage capacity up to 8 TB is furnished in the live server.

## Gallery

